

Western Blot Protocol

Required Solutions

10x PBS

80 g NaCl
2 g KCl
14.4 g Na₂HPO₄
2.4 g KH₂PO₄
dissolve into 800 ml ddH₂O
Adjust pH to 7.4
qs to 1000 ml
autoclave
make a 1x working solution

Buffer C (4°C)

20 mM Hepes-KOH pH 7.9 (from 0.5 M stock)
0.42 M KCl (from 3 M stock)
25% glycerol (from 80% stock)
0.1 mM EDTA (from 0.5 M stock pH 8.0)
5 mM MgCl₂ (from 1 M stock)
0.2% NP40 (from 10% stock)
Use autoclaved stocks for all ingredients.
Store at 4°C

Add right before use to 10 ml:

10 µl 1 M DTT (store aliquots at -20°C, thaw once)
50 µl 100 mM PMSF (in iso-propanol -20°C)
1 µl 10 mg/ml Leupeptin (-20°C)
1 µl 10 mg/ml Aprotinin (-20°C)
10 µl 1 mg/ml Pepstatin (-20°C)

Or, replace the last three protease inhibitors with one mini Roche protease inhibitor tablet (# 11 836 153001) and add DTT and PMSF as above.

For looking at phosphorylated proteins, add the following phosphatase inhibitors to Buffer C:

Na-beta-glycerophosphate to 50 mM (stock 1 M)
NaF to 1 mM (stock 0.4 M)
Na-ortho-vanadate to 1 mM (stock 0.1 M)

Buffer D

20 mM Hepes-KOH pH 7.9
100 mM KCl
25% glycerol
0.1 mM EDTA (cold)

Add right before use, for 250 ml:

250 μ l 1 M DTT

1.25 ml 100 mM PMSF

10x Ponceau S

Ponceau S 2% (w/v)

TCA 30% (w/v)

Dilute to 1x before use (can be re-used many times).

4x Laemmli Buffer

4.4 ml 0.5 M Tris (pH 6.8)

4.4 ml Glycerol

2.2 ml 20% SDS

0.5 ml 1% Bromophenol Blue

0.5 ml Beta-ME

Aliquot and store at -20°C. Dilute to 2x before use.

10x Running Buffer

30.3 g (0.25 M) Tris Base

144 g (1.92 M) Glycine

10 g (1%) SDS or appropriate for concentrated stock

qs 1000 ml ddH₂O

Dilute 1:10 with ddH₂O. pH will be 8.3

10x Towbin's Electrotransfer Buffer

30.3 g Tris Base

144 g Glycine

qs 1000 ml ddH₂O

1x Towbin's

100 ml 10x stock

200 ml (20%) Methanol

qs 1000 ml ddH₂O

1 ml (0.02%) 20% SDS— optional

Laemmli Buffer Cell Extracts

Cells can be directly lysed into Laemmli buffer as follows:

1. Harvest cells by trypsination and count cells. Keep everything cold after this step.
2. Spin cells in media for 5 min at 1,000 g at 4°C.
3. Aspirate off media and resuspend pellet in 1 ml cold 1x PBS.
4. Transfer 10^6 cells to an eppendorf tube and spin in microfuge for 4 min at 4,000 g at 4°C. Aspirate off PBS.
5. Suspend cells in 100 µl hot Laemmli buffer.
6. Heat for 5 min (80- 100°C).
7. Shear DNA through a 28 1/2 gauge insulin syringe 5 times, or sonicate.
8. Heat to 80-100°C for 5 minutes before loading.
9. Load 10 µl (10^5 cells) per lane.

Cells can also be scraped into Laemmli buffer with a rubber policeman but this procedure does not allow for a cell count.

Buffer C Whole Cell Extracts

1. Trypsinize cells (minimally one 6 well dish well or one 10 cm plate; keep cells subconfluent). Keep everything cold after this step; chill tubes and solutions on ice.
2. Wash once with medium containing serum (to inactivate the trypsin), 2 times with cold 1x PBS. Cells can be counted in one of these steps.
3. Resuspend cell pellet in 5x the pellet volume of buffer C (4°C).
4. Incubate on ice for 30 min with occasional mixing (flicking the tube).
5. Spin 15K RPM (microfuge at max setting), 10 min at 4°C.
6. Set some supernatant aside for a Bradford assay and quick freeze the sample or aliquots on dry ice and store at -80°C.
7. Optional: Dialyze to 50 vol buffer D for 2 hours at 4°C before freezing but in this case, spin dialyzed sample for 5 min in microfuge at 4°C and transfer the supernatant to a new tube before freezing.

SDS-PAGE Gel and Western Blot

Run the gel according to the instructions below. Refer to [SDS-PAGE Gels](#) on page 7 for details.

1. Measure protein concentration in duplicate using Bradford using a BSA standard curve. Run up to 70 µg/lane (expect 300-400 µg protein from a subconfluent 10 cm plate).
2. Mix extract 1:1 with Laemmli buffer and heat to 80-100°C for 5 min (don't boil too long, proteins get destroyed).
3. Run on an SDS-PAGE minigel until the blue front is at the bottom of the gel (refer to [SDS-PAGE Gels](#) for gel composition).

Bio-Rad Mini-Gel Box Running Conditions: 75 V x 3 hrs (dye about in the bottom of the gel)

4. Transfer: Blot onto a nitrocellulose membrane. Pre-wet materials in transfer buffer. Stack in the following order:

- case (clear side)
- sponge
- Whatman paper
- membrane
- gel
- Whatman paper
- sponge
- case (black side)

Place in the transfer apparatus with black side facing black).

Bio-Rad Mini-Gel Box Electrotransfer:

- 90 V x 1 hr (use ice-pack to cool down the apparatus)
- 70 V x 2 hr (use ice-pack to cool down the apparatus)

Transfer at 0.5 A-Hour (0.5 A for one hour, or 0.05 A for 10 hrs.)

For POT1 and TPP1 refer to the [Detection of POT1 and TPP1 \(mouse and human\)](#) protocol on page 5.

5. Stain with 1x Ponceau S for a minute and destain in acidified H₂O (2.3 l of ddH₂O + 4 ml concentrated HCl). Wrap in plastic wrap and Xerox. Rinse in 1x PBS.
6. Block the membrane for 30 min in 20-30 ml 1x PBS + 5% non-fat dry milk + 0.1% Tween 20, in a small Tupperware dish on a shaker.
7. Incubate with primary antibody (see [Antibodies for Shelterin Components](#) on page 8) diluted in 2 ml 1x PBS + 5% milk + 0.1% Tween 20. Incubate o/n at 4°C or 4 hrs in a seal-a-meal bag on a nutator at RT. Membranes can be stacked back to back in one bag. The primary antibody mix can be re-used (store at 4°C.)
8. Wash 3 x for 5-10 min in ~50 ml 1x PBS + 0.1% Tween 20 at RT in a small Tupperware on a shaker.
9. Incubate with secondary antibody for 30 min to 1 hr at RT in 2 ml 1x PBS + 5% milk + 0.1% Tween 20 in a seal-a-meal bag on a nutator. Use Amersham anti-mouse or rabbit HRPO conjugated at a dilution of 1:2500.
10. Wash 3 x 10 min each in ~50 ml 1x PBS + 0.1% Tween 20 at RT in a small Tupperware on a shaker.
11. Rinse with ddH₂O.
12. Detect protein with ECL kit (2 ml/membrane). In a separate tube, mix black and white ECL solutions in a 1:1 ratio. Then aliquot solution onto membranes and wait for 1 minute. Drain the ECL, wrap in plastic and expose to film.

Notes:

For ATM and p19ARF, protein must be blotted onto PVDF membranes.

For RB, proteins must be run on a 5% gel & the 50 kD MW marker should run off gel.

Stripping/Re-probing Western Blot

1. Shake filters in 2M Glycine (pH 2.2) for 20-30 minutes.
2. Rinse the filters in 1x PBS + 0.1% Tween 20 (several changes).
3. Re-block with 1x PBS + 10% non-fat dry milk + 0.1% Tween 20. The filter can now be used for the next antibody.

Note:

Some strong signals cannot be stripped of. In such case, use the membrane for probing proteins at very different molecular weight. You can strip the membrane 2-6 times.

Detection of POT1 and TPP1 (mouse and human)

Immediately after transfer:

1. Wash blot for 30 min at RT in 6 M AC buffer.
2. Wash 30 min at RT in 3 M AC buffer.
3. Wash 30 min at 4°C in 1M AC buffer.
4. Wash 30 min at 4°C 0.1 M AC buffer.
5. Wash 2 hours or overnight in AC buffer at 4°C.
6. Proceed with Western protocol, step 6 (blocking) on page 4.

AC Buffers:

(Make 8 M stock of Guanidine.)

Guanidine-HCl	6 M	3 M	1 M	0.1 M	AC
Glycerol (50%)	5 ml	5 ml	5 ml	5 ml	5 ml
5 M NaCl	0.5 ml	0.5 ml	0.5 ml	0.5 ml	0.5 ml
1 M Tris 7.5	0.5 ml	0.5 ml	0.5 ml	0.5 ml	0.5 ml
0.5 M EDTA	50 µl	50 µl	50 µl	50 µl	50 µl
Tween 20 (10%)	250 µl	250 µl	250 µl	250 µl	250 µl
8 M Guanidine-HCl	18.75 ml	9.3 ml	3.1 ml	0.31 ml	---
Milk Powder	0.5 g	0.5 g	0.5 g	0.5 g	0.5 g
1 M DTT	25 µl	25 µl	25 µl	25 µl	25 µl
ddH ₂ O	-	10.3 ml	15.6 ml	18.4 ml	18.7 ml
Total Volume:	25 ml	25 ml	25 ml	25 ml	25 ml

Note: The milk in the 1 M AC buffer does not go into solution well.

mTRF2 1254 Western Blot

The key to getting good resolution of the TRF2 band is to use a homemade gel and to run the gel slowly (65 V) until the 37 kD marker is at the bottom of the gel.

Blocking Buffer

0.1% Tween 20


5% powdered milk in 1x PBS

Incubation Buffer

0.1% Tween 20

0.1% milk in 1x PBS

1. Run gel at 65 V until the 37 kD marker is at the bottom of the gel (~3 hours).
2. Transfer the gel to nitrocellulose at 90 V for 1.5 hours (use an ice pack to cool down the apparatus) or at 100 mA overnight in the cold room.
3. (Optional) Stain blot with Ponceau to visualize bands; destain with water; make a copy on the copier (wrap in plastic wrap).
4. Block blot for 30 min at RT (in a small tupperware on a shaker).
5. Incubate blot with α mTRF2 1254 at a concentration of 1:5,000 for 2 hours at RT or overnight in the cold room on a nutator (in a seal-a-meal bag).
6. Wash blot 3 x 10 min in incubation buffer or PBST (in a small Tupperware on a shaker).
7. Incubate blot with secondary antibody (rabbit) at a concentration of 1:2,500 in incubation buffer for 30-45 min at RT on a nutator (in a seal-a-meal bag).
8. Wash the blot 3 x 10 min in incubation buffer or PBST, wash for 5 min in 1x PBS, and rinse with water (in a small Tupperware on a shaker).
9. Place blot on a piece of plastic wrap on your bench. Add ECL on top of the blot for 1 min (use equal parts of solution 1 and solution 2).
10. Shake off excess solution from the blot and wrap securely in a new piece of plastic wrap.
11. Expose blot for 1-5 min.



SDS-PAGE Gels

6% Gel	10 ml	20 ml	30 ml	60 ml
40% Acrylamide/Bis (29:1)	1.5 ml	3 ml	4.5 ml	9 ml
1.5 M Tris (pH 8.8)	2.5 ml	5 ml	7.5 ml	15 ml
ddH ₂ O	5.8 ml	11.6 ml	17.4 ml	34.8 ml
10% SDS	100 µl	200 µl	300 µl	600 µl
10% APS	100 µl	200 µl	300 µl	600 µl
TEMED	10 µl	20 µl	30 µl	60 µl

8% Gel	10 ml	20 ml	30 ml	60 ml
40% Acrylamide/Bis (29:1)	2 ml	4 ml	6 ml	12 ml
1.5 M Tris (pH 8.8)	2.5 ml	5 ml	7.5 ml	15 ml
ddH ₂ O	5.3 ml	10.6 ml	15.9 ml	31.8 ml
10% SDS	100 µl	200 µl	300 µl	600 µl
10% APS	100 µl	200 µl	300 µl	600 µl
TEMED	10 µl	20 µl	30 µl	60 µl


10% Gel	10 ml	20 ml	30 ml	60 ml
40% Acrylamide/Bis (29:1)	2.5 ml	5 ml	7.5 ml	15 ml
1.5 M Tris (pH 8.8)	2.5 ml	5 ml	7.5 ml	15 ml
ddH ₂ O	4.8 ml	9.6 ml	14.4 ml	28.8 ml
10% SDS	100 µl	200 µl	300 µl	600 µl
10% APS	100 µl	200 µl	300 µl	600 µl
TEMED	10 µl	20 µl	30 µl	60 µl

12% Gel	10 ml	20 ml	30 ml	60 ml
40% Acrylamide/Bis (29:1)	3 ml	6 ml	9 ml	18 ml
1.5 M Tris (pH 8.8)	2.5 ml	5 ml	7.5 ml	15 ml
ddH ₂ O	4.3 ml	8.6 ml	12.9 ml	25.8 ml
10% SDS	100 µl	200 µl	300 µl	600 µl
10% APS	100 µl	200 µl	300 µl	600 µl
TEMED	10 µl	20 µl	30 µl	60 µl

Stacking Gel (4%)	5 ml	10 ml	15 ml
40% Acrylamide/Bis (29:1)	0.5 ml	1.0 ml	1.5 ml
0.5 M Tris (pH 6.8) (4x)	1.25 ml	2.5 ml	3.75 ml
ddH ₂ O	3.2 ml	6.4 ml	9.6 ml
10% SDS	50 µl	100 µl	150 µl
10% APS	50 µl	100 µl	150 µl
TEMED	5 µl	10 µl	15 µl


Antibodies for Shelterin Components

Protein:	h Trf1	HeLa WCE
Antigen:	Baculovirus full length TRF1	370/371
Serum:	Rabbit Polyclonal Aff. Pur	
Number:	370/371	
Concentration:	1:1000-2000	
Reference:	Van Steensel and de Lange, 1997	



- non-specific (~70 kD)
= TRF1
-50 kD

Protein:	h Trf2	Luc TRF2 sh
Antigen:	Baculovirus full length hTRF2	
Serum:	Rabbit Polyclonal Aff. Pur	
Number:	647	
Concentration:	1:4000	
Reference:	Zhu et al, 2000	



150-
100-
75-
50-

*
*
= TRF2
*

* Non specific bands

105 HTC75 cells/lane
anti hTrf2 (647) 1:4000 1 hr RT

Protein:	h Trf2	MW
Antigen:	aa16-42 peptide	(kD)
Serum:	Rabbit Polyclonal Aff. Pur	
Number:	508	
Concentration:	1:300	
Comments:	Basic domain	
Reference:	Van Steensel et al., Cell 1998	





uninduced


0 1 2 4 6 9 days

67 -
45 -


TRF2
(clone B27)
#508




Protein: h Tin2
Antigen: Baculovirus full length hTin2
Serum: Rabbit Polyclonal Aff. Pur
Number: 864
Concentration: 1:2000
Comments: run long on a 10% gel





Protein: Tankyrase
Antigen: aa973-1149 peptide
Serum: Rabbit Polyclonal Aff. Pur
Number: 465
Concentration: 1:500
Cross reacts with: mouse, rat, hamster




Protein: m Rap1
Antigen: GST - mRap1
Serum: Rabbit Polyclonal Aff. Pur
Number: 1252
Concentration: 1:10000





Protein: m Pot1b
Antigen: aa285-307 peptide
Serum: Rabbit Polyclonal Aff. Pur
Number: 1223
Concentration: 1:1000
Comments: Guanidine renat



Protein: H Rif1
Antigen: Aa2076-2101 peptide
Serum: Rabbit Polyclonal Aff. Pur
Number: 1060
Concentration: 1:1000
Comments: Guanidine renat
Reference: Silverman et al., 2004

